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Patient Information Leaflet Revamil Melginate

Revamil Melginate 5x5 cm

The Product
Revamil Melginate – Alginate dressing met honing
5x5 cm / 10x10 cm honey-impregnated alginate sheet for the treatment of wounds

Product description

Revamil Melginate is a calcium alginate sheet impregnated with pure honey that has been produced under controlled conditions. The carefully controlled production process guarantees that Revamil Melginate with honey does not contain any traces

of pesticides and that the same quality can be reproduced consistently. When Revamil Melginate is applied on exuding wounds a gel is formed by absorption of wound fluid which protects the wound against maceration. Due to the properties of the alginate applied, i.e. the alginate is fast gelling, non-dispersible and has a high wet integrity, the alginate gel covers the wound and maintains a moist wound environment. The alginate sheet impregnated with honey, before and after gelling, protects the wound against bacterial ingress. The gel facilitates easy dressing change because it does not stick to wound and no residues remain in the wound after removal.

Instructions for use

- When required first clean the wound in the conventional way
- Open the outer packaging of the pouch
- Wear disposable gloves when handling the alginate sheet
- Remove the dressing from the transparent blister
- Apply dressing to wound surface
- Cover and secure Revamil Melginate dressing with an appropriate non-occlusive secondary dressing
- In case of high exudate level change the dressing every 24 hours
- In case of moderate exudate level change the dressing every 3 days or earlier when the dressing has reached its absorbance capacity
- If part of the alginate dressing adheres to the wound, moisten the dressing with sterile water or normal saline solution
- Gently discard the alginate dressing from the wound and discard

Indications

Revamil Melginate impregnated with honey is intended for the management of different types of exuding wounds, including bed sores, different kinds of ulcers, infected wounds, surgery and radiation induced oncological wounds, minor (1st- and 2nd degree) burns.

Contra indications

- Do not use the product if the patient is known to be allergic or hypersensitive to honey. An allergic reaction or hypersensitivity can cause redness, irritation and possibly minor pain in the skin or wound.
- Do not use on a dry wound. Using Revamil Melginate on a dry wound could lead to lower performance of the product and slower wound healing.
- In the event of unexpected side effects after use, contact your doctor or the manufacturer. Side effects can cause inconvenience to the patient.

Side effect

During the first minutes after application, the product can cause a light stinging sensation. The stinging sensation is probably caused by the low pH and high osmolality of the honey and has no effect on the performance of the product.

Safety precautions

- Do not use if honey is leaking from the inner packaging. Leakage can lead to contamination of the product with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not use if the packaging (pouch) is not tightly sealed. Packaging that is not tightly sealed can lead to non-sterility of the wound dressing and to wound infection. Because of the high antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not apply this product for more than 30 days continuously. Treatment for more than 30 days continuously could lead to inconvenience to the patient, such as irritation and/or minor pain.
- Do not reuse an already used dressing. Reuse of a dressing can lead to loss of function of the honey, contamination of the product and/or contamination or infection of the wound.
- Do not resterilize an opened packaging. Resterilization of an opened packaging can lead to loss of function of the honey, contamination of the product and/or contamination or infection of the wound.
- Do not use the product after the expiration date as this can lead to lower performance of the product and to inconvenience for the patient.

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DK: Axel Madsen A/S, Vedbæk