



**OSWELL
PENDA**
PHARMACEUTICALS

01691 674046 info@opp-uk.com

Revamil advanced wound care products embody the wound healing capacities of honey.



Revamil Melginate Wound Dressing Sizes: 5x5cm, 10x10cm

- High moisture absorption capacity for use on exuding wounds
- Enzyme rich Revamil honey offers protection against wound infections
- Melginate's fast acting gel formation helps maintain a moist wound environment
- Melginate borders are not impregnated with honey to prevent maceration
- Easy dressing change, does not adhere to the wound
- No residues remain in the wound after removal

The Revamil brand is registered for a broad spectrum of wounds.
Oswell Penda Pharmaceuticals are the UK distributor for the Revamil product range.

opp-uk.com

Revamil Melginate is a calcium alginate sheet impregnated with pure medical grade honey. The carefully controlled production process guarantees that Revamil Melginate does not contain any traces of pesticides or parabens and that the same quality can be reproduced consistently.



When Revamil Melginate is applied on exuding wounds a gel is formed by absorption of wound fluid which protects the wound against maceration.

Order Code:

Melginate	Dressings per pack	Product PIP Code
5cm x 5cm	5	<i>Code available shortly</i>
10cm x 10cm	5	<i>Code available shortly</i>

REVAMIL® PRODUCT CATEGORIES AT A GLANCE

TYPE OF WOUND	REVAMIL® GEL	REVAMIL® WOUND DRESSING	REVAMIL® BALM	REVAMIL® MELGINATE
Deep, stagnating wounds	● ●	—	—	● ● ●
Exuding wounds	—	—	—	● ● ●
Superficial wounds	●	● ● ● ●	● ● ●	—
Infected, necrotic wounds	● ● ●	—	—	● ● ●
Why, which product?	Moist wound environment Antibacterial protection	Long lasting function Hygienic Does not stick to the wound	Improves skin condition around wound For final stage of wound healing	Wound moisture management Prevents maceration